

PCR agarose

A KEY REAGENT FOR ELECTROPHORETIC SEPARATION OF DNA OR RNA FRAGMENTS

(Cat. No. P045, P051, P052)

rev. 07/2025

Description

An integral component of PCR techniques is electrophoretic separation of DNA fragments in agarose gel. However, not every agarose is suitable for preparation of gels and electrophoretic separation. PCR agarose is the agarose with optimal properties for separation of DNA and also RNA molecules. This agarose has no RNase or DNase activity and is suitable for electrophoretic separation and transfer of DNA and RNA to membranes (blotting). Low value of electroendosmosis (EEO) contributes to the increased electrophoretic mobility of RNA or DNA in the gel

Technical data

Storage

- At room temperature.

Packaging

- 100 g in plastic bottle with screw cap.

Quality control

- Each batch of PCR agarose is tested in RT-PCR.

Recommended concentration of PCR agarose for separation of DNA fragments of different sizes

| Fragment size (base pairs) | Concentration (% w/vol) of PCR agarose in gel* | |
|-------------------------------|--|---------------|
| | 1x TBE buffer | 1x TAE buffer |
| 1 000 - 23 000 | 0.60 | 0.50 |
| 800 - 10 000 | 0.80 | 0.70 |
| 400 - 8 000 | 1.00 | 0.85 |
| 300 - 7 000 | 1.20 | 1.00 |
| 200 - 4 000 | 1.50 | 1.25 |
| 100 - 3 000 | 2.00 | 1.75 |

*Final concentrations of PCR agarose in % (w/vol) in TBE or TAE buffer

Specification

| | |
|-------------------------|--------------|
| Solidification point: | 36°C |
| EEO (-m _r): | 0.09-0.13 |
| Gel force (1% g/cm): | > 1 200 |
| DNase a RNase activity: | undetectable |
| Sulphate: | < 0.15% |

| Cat. No. | Product name and specification | Amount |
|----------|--------------------------------|-----------|
| P045 | PCR agarose | 100 g |
| P051 | PCR agarose | 5x 100 g |
| P052 | PCR agarose | 10x 100 g |

