

PCR H₂O

SPECIALLY ADJUSTED WATER, RNASE AND DNASE FREE

(Cat. No P042, P242, P442)

rev. 01/2017

Description

PCR H₂O is ultrapure water (18 Mohm.cm, ultrafiltered), from which RNases were removed by treatment with diethylpyrocarbonate and by autoclaving.

Technical data

Storage

- Store at room temperature.

Packaging

- 15 ml in plastic bottle with screw cap (P042, P242).
- 50x 2.0 ml in plastic bottle with screw cap in plastic box (P442).

Quality control

- Each batch of PCR H₂O is tested in PCR.

Cat. No.	Product name and specification	Amount
P042	PCR H ₂ O	15 ml
P242	PCR H ₂ O	10x 15 ml
P442	PCR H ₂ O box	50x 2.0 ml

