

## PCR H<sub>2</sub>O

SPECIALLY ADJUSTED WATER, RNASE AND DNASE FREE

(Cat. No P042, P242)

### Description

PCR H<sub>2</sub>O is ultrapure water (18 Mohm.cm, ultrafiltered), from which RNases were removed by treatment with diethylpyrocarbonate and by autoclaving.

### Technical data

#### Storage

- Store at room temperature.

#### Packaging

- 15 ml in plastic bottle with screw cap.

#### Quality control

- Each batch of PCR H<sub>2</sub>O is tested in PCR.

Cat. No.	Product name and specification	Amount
P042	PCR H <sub>2</sub> O	15 ml
P242	PCR H <sub>2</sub> O	10x 15 ml

