

CERTIFICATE OF ANALYSIS

Product: qPCR 2x Blue Master Mix_Blue

Catalog No: B621, B622, B623, B623xl

Lot No: B621122026

Date of Expiry: 12/2026

Composition: 2x concentrated LA Hot Start Plain Master Mix contains: 150 mM Tris-HCl, pH 8.8 (at 25°C), 40 mM (NH₄)₂SO₄, 0.02% Tween 20, 5 mM MgCl₂, 400 μM dATP, 400 μM dCTP, 400 μM dGTP, 400 μM dTTP, 100 U/ml Taq DNA polymerase blend, monoclonal antibody anti-Taq (38 nM), stabilizers and additives.

Supplied with: PCR Ultra H₂O (Cat. No. P040).

Storage temperature: For short terms (days) at 4°C ± 3°C.
For long terms at -20 ± 5°C. Material can be repeatedly defrosted.

Purity: Purity of Taq DNA polymerases is verified electrophoretically (SDS PAGE). Material is free of nucleases.

Functional Test: The lot has been tested for the ability to amplify a fragment of genomic DNA using the following conditions:

Test conditions:

Volume*	Reagent	Final concentration
12.5 μl	qPCR 2x Blue	1x qPCR 2x Blue
0.5 μl	Forward primer	50 μM 5' primer 5'-ATGAACCCAGCCATCAGCG-3'
0.5 μl	Reverse primer	50 μM 3' primer 5'-GGGTAAGGACCTTGATATAGG-3'
1 μl	Template DNA	containing 80 ng of mouse genomic DNA
10.5 μl	PCR Ultra H ₂ O	(to a final volume 25 μl)

Cycling conditions:

	Temperature	Time	Number of cycles
Initial denaturation	94°C	1 min	1
Denaturation	94°C	15 s	30
Annealing of primers	55°C	15 s	
Extension	72°C	1 min	
Final extension	72°C	7 min	1
Cooling	22°C		

Result: As expected, electrophoresis of the PCR product on agarose gel revealed one band of 864 bp

FOR RESEARCH USE**APPROVED DATE:** 10.10.2024

Manager: Hana Těšitelová