

CERTIFICATE OF ANALYSIS

Product:	Taq DNA polymeráza 1.1
Catalog No:	T112, T113, T114
Lot No:	T112122027
Date of Expiry:	12/2027
Concentration:	1U/μl
Storage buffer:	20 mM Tris-HCl (pH 8.0 at 25°C), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Nonidet P-40, 0.5% Tween 20, inert red dye, stabilizers, 50% glycerol.
Supplied with:	10x PCR Blue Buffer: 750 mM Tris-HCl, pH 8.8 (at 25°C), 200 mM (NH ₄) ₂ SO ₄ , 1% Tween 20, 25 mM MgCl ₂ .
Storage temperature:	-16 to -25 °C
Purity:	The enzyme was analyzed by SDS-PAGE and single band of ~94 kDa was observed
Functional Test:	The Lot has been tested for the ability to amplify a fragment of genomic DNA using the following conditions: Test conditions: 39.5 μl PCR H ₂ O 5 μl 10x PCR Blue Buffer with 25 mM MgCl ₂ 1 μl 10 mM dNTP mix (10 mM for each, dATP, dCTP, dGTP, and d TTP) 0.5 μl 50 μM 5' primer (5'-ATGAACCCAGCCATCAGCG-3') 0.5 μl 50 μM 3' primer (5'-GGGTAAGGACCTTGATATAGG-3') 2.5 μl Taq DNA polymeráza 1.1 (2.5 U total) 1 μl DNA containing 80 ng of mouse genomic (tail) DNA. Cycling conditions: 95°C, 2 min initial denaturation, followed by 40 cycles of 94°C, 15 s (denaturation) 54°C, 15 s (annealing) 72°C, 60 s (extension) Result: As expected, electrophoresis of the PCR product on agarose gel revealed one band of 864 bp

FOR RESEARCH USE

APPROVED DATE: 12.09.2025

Manager: Hana Těšitelová