

PCR Ethanol, 75%75% ETHANOL IN PCR H₂O

(Cat. No. P044)

rev. 04/2025

Description

PCR Ethanol, 75% is 75% ethyl alcohol in ultrapure water (18 Mohm.cm, ultrafiltered), from which RNases were removed by the treatment with diethylpyrocarbonate followed by autoclaving.

Technical data**Storage**

- Store at temperature 4°C ± 3°C.

Packaging

- 100 ml in glass bottle with screw cap.

Quality control

- Each batch of PCR Ethanol, 75% is tested in RT-PCR.

Cat. No.	Product name and specification	Amount
P044	PCR Ethanol, 75%	100 ml

